

兔抗 SDC4(Ab-179) 多克隆抗体

中文名称：兔抗 SDC4(Ab-179) 多克隆抗体

英文名称： Anti-SDC4(Ab-179) rabbit polyclonal antibody

别名： syndecan 4

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SDC4(Ab-179)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

The protein encoded by this gene is a transmembrane (type I) heparan sulfate proteoglycan that functions as a receptor in intracellular signaling. The encoded protein is found as a homodimer and

	is a member of the syndecan proteoglycan family. This gene is found on chromosome 20, while a pseudogene has been found on chromosome 22.
Applications:	WB, IHC, IF
Name of antibody:	SDC4(Ab-179)
Immunogen:	Synthesized non-phosphopeptide derived from human Syndecan4 around the phosphorylation site of serine 179 (E-G-S(p)-Y-D).
Full name:	syndecan 4
Synonyms :	syndecan 4
SwissProt:	P31431
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	22 kDa
WB Positive control:	HepG2 cells lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells

IF Recommend dilution:

100-500



