

## 兔抗 PDCD5 多克隆抗体

中文名称：兔抗 PDCD5 多克隆抗体

英文名称： Anti-PDCD5 rabbit polyclonal antibody

别名： programmed cell death 5; TFAR19

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDCD5

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is upregulated during apoptosis where it translocates rapidly from the cytoplasm to the nucleus. The encoded protein may be an important regulator of K(lysine) acetyltransferase 5 (a protein involved in transcription, DNA damage response and cell cycle control) by inhibiting its proteasome-dependent degradation. Pseudogenes have been identified on chromosomes 5 and 12
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PDCD5
<b>Immunogen:</b>	Fusion protein of human PDCD5

<b>Full name:</b>	programmed cell death 5
<b>Synonyms:</b>	TFAR19
<b>SwissProt:</b>	O14737
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-200

