

兔抗 SDCCAG8 多克隆抗体

中文名称：兔抗 SDCCAG8 多克隆抗体

英文名称： Anti-SDCCAG8 rabbit polyclonal antibody

别名： BBS16; CCCAP; SLSN7; NPHP10; hCCCAP; HSPC085; NY-CO-8; CCCAP SLSN7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SDCCAG8

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a centrosome associated protein. This protein may be involved in organizing the centrosome during interphase and mitosis. Mutations in this gene are associated with retinal-renal ciliopathy.
Applications:	ELISA, IHC
Name of antibody:	SDCCAG8
Immunogen:	Synthetic peptide of human SDCCAG8
Full name:	serologically defined colon cancer antigen 8
Synonyms :	BBS16; CCCAP; SLSN7; NPHP10; hCCCAP; HSPC085 ; NY-CO-8; CCCAP SLSN7

SwissProt:	Q86SQ7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and human lung cancer
IHC Recommend dilution:	25-100



