

## 兔抗 PDCD6IP 多克隆抗体

中文名称：兔抗 PDCD6IP 多克隆抗体

英文名称： Anti-PDCD6IP rabbit polyclonal antibody

别名： AIP1; ALIX; HP95; DRIP4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDCD6IP

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a protein that functions within the ESCRT pathway in the abscission stage of cytokinesis, in intraluminal endosomal vesicle formation, and in enveloped virus budding. Studies using mouse cells have shown that overexpression of this protein can block apoptosis. In addition, the product of this gene binds to the product of the PDCD6 gene, a protein required for apoptosis, in a calcium-dependent manner. This gene product also binds to endophilins, proteins that regulate membrane shape during endocytosis. Overexpression of this gene product and endophilins results in cytoplasmic vacuolization, which

	may be partly responsible for the protection against cell death. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. Related pseudogenes have been identified on chromosome 15.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	PDCD6IP
<b>Immunogen:</b>	Synthetic peptide of human PDCD6IP
<b>Full name:</b>	programmed cell death 6 interacting protein
<b>Synonyms:</b>	AIP1; ALIX; HP95; DRIP4
<b>SwissProt:</b>	Q8WUM4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	96 kDa
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	200-1000

