

## 兔抗 PDCD6 多克隆抗体

中文名称：兔抗 PDCD6 多克隆抗体

英文名称： Anti-PDCD6 rabbit polyclonal antibody

别名： ALG-2, PEF1B

抗原： PDCD6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. Calcium binding is important for homodimerization and for conformational changes required for binding to other protein partners. This gene product participates in T cell receptor-, Fas-, and glucocorticoid-induced programmed cell death. In mice deficient for this gene product, however, apoptosis was not blocked suggesting this gene product is functionally redundant. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is also located on the short arm

	of chromosome 5.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PDCD6
<b>Immunogen:</b>	Fusion protein of human PDCD6
<b>Full name:</b>	programmed cell death 6
<b>Synonyms :</b>	ALG-2, PEF1B
<b>SwissProt:</b>	O75340
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human ovarian cancer and human colon cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Jurkat and hela cells
<b>WB Recommended dilution:</b>	200-1000



