

兔抗 PDCL2 多克隆抗体

中文名称：兔抗 PDCL2 多克隆抗体

英文名称：Anti-PDCL2 rabbit polyclonal antibody

别名：phosducin like 2; GCPHLP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDCL2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the phosducin-like protein family and is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin. Members of the phosducin-like protein family have been shown to bind to the beta-gamma subunits of G proteins.
Applications:	ELISA, IHC
Name of antibody:	PDCL2
Immunogen:	Fusion protein of human PDCL2
Full name:	phosducin like 2
Synonyms:	GCPHLP

SwissProt:	Q8N4E4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300

