

兔抗 VPS37A 多克隆抗体

中文名称: 兔抗 VPS37A 多克隆抗体

英文名称: Anti-VPS37A rabbit polyclonal antibody

别名: HCRP1; PQBP2; SPG53

储存: 冷冻 (-20℃)

抗原: VPS37A

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene belongs to the VPS37 family, and encodes a component of the ESCRT-I (endosomal sorting complex required for transport I) protein complex, required for the sorting of ubiquitinated transmembrane proteins into internal vesicles of multivesicular bodies. Expression of this gene is downregulated in hepatocellular carcinoma, and mutations in this gene are associated with autosomal recessive spastic paraplegia-53. A related pseudogene has been identified on chromosome 5. Alternatively

	spliced transcript variants have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	VPS37A
Immunogen:	Fusion protein of human VPS37A
Full name:	VPS37A subunit of ESCRT-I
Synonyms:	HCRP1; PQBP2; SPG53
SwissProt:	Q8NEZ2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	44 kDa
WB Positive control:	Human bladder carcinoma and fetal liver tissue
WB Recommended dilution:	500-2000

