

## 兔抗 PDCL 多克隆抗体

中文名称：兔抗 PDCL 多克隆抗体

英文名称： Anti-PDCL rabbit polyclonal antibody

别名： phosducin like; PhLP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDCL

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Phosducin-like protein is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin, a phosphoprotein expressed in retina and pineal gland. Both phosducin-like protein and phosphoducin have been shown to regulate G-protein signaling by binding to the beta-gamma subunits of G proteins.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PDCL
<b>Immunogen:</b>	Fusion protein of human PDCL
<b>Full name:</b>	phosducin like

<b>Synonyms:</b>	PhLP
<b>SwissProt:</b>	Q13371
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-50

