

## 兔抗 VPS37D 多克隆抗体

中文名称: 兔抗 VPS37D 多克隆抗体

英文名称: Anti-VPS37D rabbit polyclonal antibody

别名: WBSCR24

抗原: VPS37D

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies. May be involved in cell growth and differentiation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	VPS37D
<b>Immunogen:</b>	Fusion protein of human VPS37D
<b>Full name:</b>	VPS37D subunit of ESCRT-I

<b>Synonyms:</b>	WBSCR24
<b>SwissProt:</b>	Q86XT2
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

