

## 兔抗 SDF4 多克隆抗体

中文名称：兔抗 SDF4 多克隆抗体

英文名称： Anti-SDF4 rabbit polyclonal antibody

别名： Cab45; SDF-4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SDF4

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a stromal cell derived factor that is a member of the CREC protein family. The encoded protein contains six EF-hand motifs and calcium-binding motifs. This protein localizes to the Golgi lumen and may be involved in regulating calcium dependent cellular activities.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SDF4
<b>Immunogen:</b>	Fusion protein of human SDF4
<b>Full name:</b>	stromal cell derived factor 4
<b>Synonyms :</b>	Cab45; SDF-4

<b>SwissProt:</b>	Q9BRK5
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human liver cancer and human prostate cancer
<b>IHC Recommend dilution:</b>	30-150

