

## 兔抗 PDE4DIP 多克隆抗体

中文名称：兔抗 PDE4DIP 多克隆抗体

英文名称： Anti-PDE4DIP rabbit polyclonal antibody

别名： phosphodiesterase 4D interacting protein; MMGL; CMYA2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDE4DIP

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene serves to anchor phosphodiesterase 4D to the Golgi/centrosome region of the cell. Defects in this gene may be a cause of myeloproliferative disorder (MBD) associated with eosinophilia. Several transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PDE4DIP
<b>Immunogen:</b>	Fusion protein of human PDE4DIP
<b>Full name:</b>	phosphodiesterase 4D interacting protein
<b>Synonyms:</b>	MMGL; CMYA2

<b>SwissProt:</b>	Q5VU43
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-200

