

兔抗 SDHA 多克隆抗体

中文名称：兔抗 SDHA 多克隆抗体

英文名称： Anti-SDHA rabbit polyclonal antibody

别名： FP; PGL5; SDH1; SDH2; SDHF; CMD1GG

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SDHA

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a major catalytic subunit of succinate-ubiquinone oxidoreductase, a complex of the mitochondrial respiratory chain. The complex is composed of four nuclear-encoded subunits and is localized in the mitochondrial inner membrane. Mutations in this gene have been associated with a form of mitochondrial respiratory chain deficiency known as Leigh Syndrome. A pseudogene has been identified on chromosome 3q29.
Applications:	ELISA, WB, IHC
Name of antibody:	SDHA

Immunogen:	Fusion protein of human SDHA
Full name:	succinate dehydrogenase complex, subunit A, flavoprotein (Fp)
Synonyms :	FP; PGL5; SDH1; SDH2; SDHF; CMD1GG
SwissProt:	P31040
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and Human thyroid cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	73 kDa
WB Positive control:	Mouse brain tissue and Hela cells, Jurkat cells and Mouse heart tissue
WB Recommended dilution:	500-2000



