

## 兔抗 VPS53 多克隆抗体

中文名称: 兔抗 VPS53 多克隆抗体

英文名称: Anti-VPS53 rabbit polyclonal antibody

别名: HCCS1; PCH2E; hVps53L; pp13624

储存: 冷冻 (-20℃)

抗原: VPS53

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with sequence similarity to the yeast Vps53p protein. Vps53p is involved in retrograde vesicle trafficking in late Golgi.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	VPS53
<b>Immunogen:</b>	Fusion protein of human VPS53
<b>Full name:</b>	VPS53 subunit of GARP complex
<b>Synonyms:</b>	HCCS1; PCH2E; hVps53L; pp13624
<b>SwissProt:</b>	Q5VIR6
<b>ELISA Recommended dilution:</b>	1000-2000

<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	80 kDa
<b>WB Positive control:</b>	A172 cell
<b>WB Recommended dilution:</b>	200-1000



