

兔抗 PDE6D 多克隆抗体

中文名称：兔抗 PDE6D 多克隆抗体

英文名称： Anti-PDE6D rabbit polyclonal antibody

别名： phosphodiesterase 6D; PDED; JBTS22

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDE6D

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the delta subunit of rod-specific photoreceptor phosphodiesterase (PDE), a key enzyme in the phototransduction cascade. A similar protein in cow functions in solubilizing membrane-bound PDE. In addition to its role in the PDE complex, the encoded protein is thought to bind to prenyl groups of proteins to target them to subcellular organelles called cilia. Mutations in this gene are associated with Joubert syndrome-22. Alternative splicing results in multiple splice variants.
Applications:	ELISA, IHC
Name of antibody:	PDE6D

Immunogen:	Fusion protein of human PDE6D
Full name:	phosphodiesterase 6D
Synonyms:	PDED; JBTS22
SwissProt:	O43924
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-100

