

## 兔抗 PDGFA 多克隆抗体

- 中文名称： 兔抗 PDGFA 多克隆抗体
- 英文名称： Anti-PDGFA rabbit polyclonal antibody
- 别名： PDGF1, PDGF-A
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 抗原： PDGFA
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene is a member of the platelet-derived growth factor family. The four members of this family are mitogenic factors for cells of mesenchymal origin and are characterized by a motif of eight cysteines. This gene product can exist either as a homodimer or as a heterodimer with the platelet-derived growth factor beta polypeptide, where the dimers are connected by disulfide bonds. Studies using knockout mice have shown cellular defects in oligodendrocytes, alveolar smooth muscle cells, and Leydig cells in the testis; knockout mice die

	either as embryos or shortly after birth. Two splice variants have been identified for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PDGFA
<b>Immunogen:</b>	Synthetic peptide of human PDGFA
<b>Full name:</b>	Platelet-derived growth factor alpha polypeptide
<b>Synonyms :</b>	PDGF1, PDGF-A
<b>SwissProt:</b>	P04085
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human colon cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	100-300

