

兔抗 PDGFD 多克隆抗体

中文名称：兔抗 PDGFD 多克隆抗体

英文名称：Anti-PDGFD rabbit polyclonal antibody

别名：platelet derived growth factor D; IEGF; SCDGFB; MSTP036; SCDGF-B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PDGFD

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the platelet-derived growth factor family. The four members of this family are mitogenic factors for cells of mesenchymal origin and are characterized by a core motif of eight cysteines, seven of which are found in this factor. This gene product only forms homodimers and, therefore, does not dimerize with the other three family members. It differs from alpha and beta members of this family in having an unusual N-terminal domain, the CUB domain. Two splice variants have been identified for this gene. [provided by RefSeq, Jul 2008]

Applications:	ELISA, IHC
Name of antibody:	PDGFD
Immunogen:	Fusion protein of human PDGFD
Full name:	platelet derived growth factor D
Synonyms:	IEGF; SCDGFB; MSTP036; SCDGF-B
SwissProt:	Q9GZP0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	100-200



