

兔抗 ITIH4 多克隆抗体

中文名称：兔抗 ITIH4 多克隆抗体

英文名称：Anti-ITIH4 rabbit polyclonal antibody

别名：inter-alpha-trypsin inhibitor heavy chain family member 4; H4P; IHRP; GP120; PK120; ITIHL1; PK-120; ITI-HC4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ITIH4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is secreted into the blood, where it is cleaved by plasma kallikrein into two smaller forms. Expression of this gene has been detected only in liver, and it seems to be upregulated during surgical trauma. This gene is part of a cluster of similar genes on chromosome 3. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ITIH4

Immunogen:	Synthetic peptide of human ITIH4
Full name:	inter-alpha-trypsin inhibitor heavy chain family member 4
Synonyms:	H4P; IHRP; GP120; PK120; ITIHL1; PK-120; ITI-HC4
SwissProt:	Q14624
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution :	30-150

