

兔抗 SELENBP1 多克隆抗体

中文名称：兔抗 SELENBP1 多克隆抗体

英文名称： Anti-SELENBP1 rabbit polyclonal antibody

别名： selenium binding protein 1; LPSB; SP56; hSBP; SBP56; HEL-S-134P

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SELENBP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the selenium-binding protein family. Selenium is an essential nutrient that exhibits potent anticarcinogenic properties, and deficiency of selenium may cause certain neurologic diseases. The effects of selenium in preventing cancer and neurologic diseases may be mediated by selenium-binding proteins, and decreased expression of this gene may be associated with several types of cancer. The encoded protein may play a selenium-dependent role in ubiquitination/deubiquitination-mediated protein degradation. Alternatively spliced transcript variants encoding multiple isoforms have been observed fo

	r this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	SELENBP1
Immunogen:	Fusion protein of human SELENBP1
Full name:	selenium binding protein 1
Synonyms:	LPSB; SP56; hSBP; SBP56; HEL-S-134P
SwissProt:	Q13228
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	53 kDa
WB Positive control:	Human fetal liver tissue and Mouse liver tissue lysates
WB Recommended dilution:	200-1000



