

兔抗 PDGFRA(Phospho-Tyr849) 多克隆抗体

中文名称：兔抗 PDGFRA(Phospho-Tyr849) 多克隆抗体

英文名称： Anti-PDGFRA(Phospho-Tyr849) rabbit polyclonal antibody

别名： CD140A; PDGFR2; PDGFR-2; RHEPDGFRA

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PDGFRA(Phospho-Tyr849)

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides.
Applications:	WB
Name of antibody:	PDGFRA(Phospho-Tyr849)
Immunogen:	Peptide sequence around phosphorylation site of Tyrosine849(S-N-Y(p)-V-S) derived from Human PDGFRA.

Full name:	platelet-derived growth factor receptor, alpha polypeptide
Synonyms :	CD140A; PDGFR2; PDGFR-2; RHEPDGFRA
SwissProt:	P16234
WB Predicted band size:	123 kDa
WB Positive control:	293 cells lysates
WB Recommended dilution:	500-1000

