

兔抗 ITK 多克隆抗体

中文名称： 兔抗 ITK 多克隆抗体

英文名称： Anti-ITK rabbit polyclonal antibody

别名： EMT; LYK; LPFS1; PSCTK2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： ITK

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation.
Applications:	ELISA, IHC
Name of antibody:	ITK
Immunogen:	Fusion protein of human ITK
Full name:	IL2-inducible T-cell kinase
Synonyms :	EMT; LYK; LPFS1; PSCTK2

SwissProt:	Q08881
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	25-100

