

兔抗 PDGFRB(Phospho-Tyr740) 多克隆抗体

中文名称：兔抗 PDGFRB(Phospho-Tyr740) 多克隆抗体

英文名称： Anti-PDGFRB(Phospho-Tyr740) rabbit polyclonal antibody

别名： IMF1; IBGC4; JTK12; PDGFR; CD140B; PDGFR1; PDGFR-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PDGFRB(Phospho-Tyr740)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides.
Applications:	WB
Name of antibody:	PDGFRB(Phospho-Tyr740)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine

	740 (G-G-Y(p)-M-D) derived from Human PDGFR beta.
Full name:	platelet-derived growth factor receptor, beta polypeptide
Synonyms :	IMF1; IBGC4; JTK12; PDGFR; CD140B; PDGFR1; PDGFR-1
SwissProt:	P09619
WB Predicted band size:	124 kDa
WB Positive control:	HepG2 cells lysates treated with EGF
WB Recommended dilution:	500-1000

