

## 兔抗 SELE 多克隆抗体

中文名称：兔抗 SELE 多克隆抗体

英文名称： Anti-SELE rabbit polyclonal antibody

别名： ELAM; ESEL; CD62E; ELAM1; LECAM2

抗原： SELE

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

The protein encoded by this gene is found in cytokine-stimulated endothelial cells and is thought to be responsible for the accumulation of blood leukocytes at sites of inflammation by mediating the adhesion of cells to the vascular lining. It exhibits structural features such as the presence of lectin- and EGF-like domains followed by short consensus repeat (SCR) domains that contain 6 conserved cysteine residues. These proteins are part of the selectin family of cell adhesion molecules. Adhesion molecules participate in the interaction between leukocytes and the endothelium and appear to be involved in the p

	athogenesis of atherosclerosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SELE
<b>Immunogen:</b>	Fusion protein of human SELE
<b>Full name:</b>	selectin E
<b>Synonyms :</b>	ELAM; ESEL; CD62E; ELAM1; LECAM2
<b>SwissProt:</b>	P16581
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200

