

兔抗 ITPKC 多克隆抗体

- 中文名称: 兔抗 ITPKC 多克隆抗体
- 英文名称: Anti-ITPKC rabbit polyclonal antibody
- 别名: inositol-trisphosphate 3-kinase C; IP3KC; IP3-3KC
- 抗原: ITPKC
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 相关类别: 一抗
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the inositol 1,4,5-trisphosphate [Ins(1,4,5)P(3)] 3-kinase family of enzymes that catalyze the phosphorylation of inositol 1,4,5-trisphosphate to 1,3,4,5-tetrakisphosphate. The encoded protein is localized to the nucleus and cytoplasm and has both nuclear import and nuclear export activity. Single nucleotide polymorphisms in this gene are associated with Kawasaki disease.
Applications:	ELISA, IHC
Name of antibody:	ITPKC
Immunogen:	Synthetic peptide of human ITPKC

Full name:	inositol-trisphosphate 3-kinase C
Synonyms:	IP3KC; IP3-3KC
SwissProt:	Q96DU7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

