

兔抗 SELL 多克隆抗体

中文名称: 兔抗 SELL 多克隆抗体

英文名称: Anti-SELL rabbit polyclonal antibody

别名: selectin L; TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; LECAM1

抗原: SELL

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a cell surface adhesion molecule that belongs to a family of adhesion/homing receptors. The encoded protein contains a C-type lectin-like domain, a calcium-binding epidermal growth factor-like domain, and two short complement-like repeats. The gene product is required for binding and subsequent rolling of leucocytes on endothelial cells, facilitating their migration into secondary lymphoid organs and inflammation sites. Single-nucleotide polymorphisms in this gene have been associated with various diseases including immunoglobulin A nephropathy. Alternatively spliced transcript variants have b

	een found for this gene.
Applications:	ELISA, IHC
Name of antibody:	SELL
Immunogen:	Synthetic peptide of human SELL
Full name:	selectin L
Synonyms:	TQ1; LAM1; LEU8; LNHR; LSEL; CD62L; LYAM1; PLNHR; L ECAM1
SwissProt:	P14151
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	20-100



