

兔抗 VSIG10 多克隆抗体

中文名称：兔抗 VSIG10 多克隆抗体

英文名称：Anti-VSIG10 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：VSIG10

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

It is the genetic home for a number of immunoglobulin superfamily members including the killer cell and leukocyte Ig-like receptors, a number of ICAMs, the CEACAM and PSG family, and Fc receptors. Key genes for eye color and hair color also map to chromosome 19. Peutz-Jeghers syndrome, spinocerebellar ataxia type 6, the stroke disorder CADASIL, hypercholesterolemia and insulin-dependent diabetes have been linked to chromosome 19. Translocations with chromosome 19 and chromosome 14 can be seen in some lymphoproliferative disorders and typically involve the proto-oncogene BCL3. The LOC147645 gene product has been provisionally designated LOC147645 pending further characteri

	zation.
Applications:	ELISA, IHC
Name of antibody:	VSIG10
Immunogen:	Synthetic peptide of humanVSIG10
Full name:	V-set and immunoglobulin domain containing 10
SwissProt:	Q8N0Z9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

