

兔抗 ITSN2 多克隆抗体

中文名称: 兔抗 ITSN2 多克隆抗体

英文名称: Anti-ITSN2 rabbit polyclonal antibody

别名: intersectin 2; SWA; SWAP; SH3D1B; SH3P18; PRO2015

抗原: ITSN2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a cytoplasmic protein which contains SH3 domains. This protein is a member of a family of proteins involved in clathrin-mediated endocytosis. Intersectin 2 is thought to regulate the formation of clathrin-coated vesicles and also may function in the induction of T cell antigen receptor (TCR) endocytosis.
Applications:	ELISA, IHC
Name of antibody:	ITSN2

Immunogen:	Synthetic peptide of human ITSN2
Full name:	intersectin 2
Synonyms:	SWA; SWAP; SH3D1B; SH3P18; PRO2015
SwissProt:	Q9NZM3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human liver cancer
IHC Recommend dilution:	40-200





