

## 兔抗 SEMA3D 多克隆抗体

中文名称：兔抗 SEMA3D 多克隆抗体

英文名称：Anti-SEMA3D rabbit polyclonal antibody

别名：coll-2; Sema-Z2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SEMA3D

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene belongs to the semaphorin family. Induces the collapse and paralysis of neuronal growth cones. Could potentially act as repulsive cues toward specific neuronal populations. Strong binding to neuropilin is mediated by the carboxy third of the protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SEMA3D
<b>Immunogen:</b>	Synthetic peptide of human SEMA3D
<b>Full name:</b>	sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3D

<b>Synonyms :</b>	coll-2; Sema-Z2
<b>SwissProt:</b>	O95025
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

