

兔抗 SEMA3G 多克隆抗体

中文名称：兔抗 SEMA3G 多克隆抗体

英文名称：Anti-SEMA3G rabbit polyclonal antibody

别名： sem2

抗原： SEMA3G

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The transcription of this gene is activated by PP AR-gamma, and the resulting protein product plays a role in endothelial cell migration. Expression of this gene also inhibits tumor cell migration and invasion.
Applications:	ELISA, WB, IHC
Name of antibody:	SEMA3G
Immunogen:	Fusion protein of human SEMA3G
Full name:	semaphorin 3G
Synonyms:	sem2
SwissProt:	Q9NS98

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human tonsil
IHC Recommend dilution:	150-300
WB Predicted band size:	87 kDa
WB Positive control:	Mouse kidney tissue lysate
WB Recommended dilution:	500-2000



