

兔抗 VSX1 多克隆抗体

中文名称：兔抗 VSX1 多克隆抗体

英文名称：Anti-VSX1 rabbit polyclonal antibody

别名：PPD; KTCN; PPCD; RINX; KTCN1; PPCD1; CAASDS

相关类别：一抗

抗原：VSX1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a paired-like homeodomain and binds to the core of the locus control region of the red/green visual pigment gene cluster. The encoded protein may regulate expression of the cone opsin genes early in development. Mutations in this gene can cause posterior polymorphous corneal dystrophy and keratoconus. Alternatively spliced transcript variants encoding different isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	VSX1

Immunogen:	Fusion protein of human VSX1
Full name:	visual system homeobox 1
Synonyms:	PPD; KTCN; PPCD; RINX; KTCN1; PPCD1; CAASDS
SwissProt:	Q9NZR4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human gastric cancer
IHC Recommend dilution:	25-100



