

兔抗 VTA1 多克隆抗体

中文名称：兔抗 VTA1 多克隆抗体

英文名称： Anti-VTA1 rabbit polyclonal antibody

别名： vesicle trafficking 1; DRG1; LIP5; SBP1; DRG-1; My012; C6orf55; HSPC228

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： VTA1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	C6ORF55 encodes a protein involved in trafficking of the multivesicular body, an endosomal compartment involved in sorting membrane proteins for degradation in lysosomes (Ward et al., 2005 [PubMed 15644320]).
Applications:	ELISA, IHC
Name of antibody:	VTA1
Immunogen:	Fusion protein of human VTA1
Full name:	vesicle trafficking 1
Synonyms:	DRG1; LIP5; SBP1; DRG-1; My012; C6orf55; HSPC228
SwissProt:	Q9NP79
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	50-200



