


## 兔抗 PDIA5 多克隆抗体

- 中文名称：兔抗 PDIA5 多克隆抗体
- 英文名称：Anti-PDIA5 rabbit polyclonal antibody
- 别名：protein disulfide isomerase family A member 5; PDIR
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：PDIA5
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, three catalytically active thioredoxin (TRX) domains, a TRX-like domain, and a C-terminal ER-retention sequence. The N-terminal TRX-like domain is the primary binding site for the major ER chaperone calreticulin and possibly other proteins and substrates as well. Alternative splicing results in multiple protein- and non-protein-coding transcript variants.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	PDIA5
<b>Immunogen:</b>	Fusion protein of human PDIA5
<b>Full name:</b>	protein disulfide isomerase family A member 5
<b>Synonyms:</b>	PDIR
<b>SwissProt:</b>	Q14554
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	60 kDa
<b>WB Positive control:</b>	Human placenta tissue lysate
<b>WB Recommended dilution:</b>	1000-5000



**mlbio** 梅联生物  
Good elisakit producers

