

兔抗 SEMA4C 多克隆抗体

中文名称：兔抗 SEMA4C 多克隆抗体

英文名称： Anti-SEMA4C rabbit polyclonal antibody

别名： SEMAF; SEMAI; SEMACL1; M-SEMA-F

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEMA4C

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Cell surface receptor for PLXNB2 that plays an important role in cell-cell signaling. PLXNB2 binding promotes downstream activation of RHOA and phosphorylation of ERB B2 at 'Tyr-1248'. Required for normal brain development, axon guidance and cell migration (By similarity). Probabl e signaling receptor which may play a role in myogenic differentiation through activation of the stress-activated MAPK cascade.
Applications:	ELISA, IHC
Name of antibody:	SEMA4C
Immunogen:	Fusion protein of human SEMA4C

Full name:	sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 4C
Synonyms :	SEMAF; SEMAI; SEMACL1; M-SEMA-F
SwissProt:	Q9C0C4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	25-100

