

兔抗 JAK2(Ab-570) 多克隆抗体

中文名称: 兔抗 JAK2(Ab-570) 多克隆抗体

英文名称: Anti-JAK2(Ab-570) rabbit polyclonal antibody

别名: JTK10; THCYT3

抗原: JAK2(Ab-570)

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

相关类别: 一抗

反应种属: Human Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.
Applications:	WB, IF
Name of antibody:	JAK2(Ab-570)
Immunogen:	Synthesized non-phosphopeptide derived from human JAK2

	around the phosphorylation site of tyrosine 570 (G-D-Y(p)-G-Q).
Full name:	Janus kinase 2
Synonyms :	JTK10; THCYT3
SwissProt:	O60674
WB Predicted band size:	131 kDa
WB Positive control:	293 cells treated with etoposide (25uM, 24hours) lysate
WB Recommended dilution:	500-3000
IF positive control:	NIH/3T3 cells
IF Recommend dilution:	100-500



