

兔抗 SEMA4F 多克隆抗体

中文名称：兔抗 SEMA4F 多克隆抗体

英文名称： Anti-SEMA4F rabbit polyclonal antibody

别名： S4F; SEMAM; SEMAW; M-SEMA; PRO2353; m-Sema-M

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEMA4F

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane class IV semaphorin family protein, which plays a role in neural development. This gene may be involved in neurogenesis in prostate cancer, the development of neurofibromas, and breast cancer tumorigenesis. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	SEMA4F
Immunogen:	Fusion protein of human SEMA4F
Full name:	sema domain, immunoglobulin domain (Ig), transmembrane domain (TM) and short cytoplasmic domain,

	(semaphorin) 4F
Synonyms :	S4F; SEMAM; SEMAW; M-SEMA; PRO2353; m-Sema-M
SwissProt:	O95754
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and human brain
IHC Recommend dilution:	25-100

