

## 兔抗 PDK2 多克隆抗体

中文名称：兔抗 PDK2 多克隆抗体

英文名称： Anti-PDK2 rabbit polyclonal antibody

别名： PDHK2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDK2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the pyruvate dehydrogenase kinase family. The encoded protein phosphorylates pyruvate dehydrogenase, down-regulating the activity of the mitochondrial pyruvate dehydrogenase complex. Overexpression of this gene may play a role in both cancer and diabetes. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	PDK2
<b>Immunogen:</b>	Synthetic peptide of human PDK2
<b>Full name:</b>	Pyruvate dehydrogenase kinase, isozyme 2

<b>Synonyms :</b>	PDHK2
<b>SwissProt:</b>	Q15119
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Mouse heart tissue
<b>WB Recommended dilution:</b>	200-1000

