

兔抗 PDLIM4 多克隆抗体

中文名称：兔抗 PDLIM4 多克隆抗体

英文名称： Anti-PDLIM4 rabbit polyclonal antibody

别名： RIL

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDLIM4

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis.
Applications:	ELISA, WB, IHC
Name of antibody:	PDLIM4
Immunogen:	Fusion protein of human PDLIM4
Full name:	PDZ and LIM domain 4
Synonyms :	RIL
SwissProt:	P50479
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer and human lung cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	35 kDa
WB Positive control:	A431 cells
WB Recommended dilution:	200-1000



