

兔抗 VTN 多克隆抗体

中文名称：兔抗 VTN 多克隆抗体

英文名：Anti-VTN rabbit polyclonal antibody

别名：VN; V75; VNT

相关类别：一抗

储存：冷冻（-20℃）

抗原：VTN

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型 rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the pexin family. It is found in serum and tissues and promotes cell adhesion and spreading, inhibits the membrane-damaging effect of the terminal cytolytic complement pathway, and binds to several serpin serine protease inhibitors. It is a secreted protein and exists in either a single chain form or a clipped, two chain form held together by a disulfide bond.
Applications:	ELISA, IHC
Name of antibody:	VTN
Immunogen:	Fusion protein of human VTN
Full name:	vitronectin
Synonyms :	VN; V75; VNT

SwissProt:	P04004
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

