

兔抗 SEMA6A 多克隆抗体

中文名称：兔抗 SEMA6A 多克隆抗体

英文名称： Anti-SEMA6A rabbit polyclonal antibody

别名： VIA; SEMA; HT018; SEMAQ; SEMA6A1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEMA6A

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The transmembrane semaphorin SEMA6A is expressed in developing neural tissue and is required for proper development of the thalamocortical projection. Cell surface receptor for PLXNA2 that plays an important role in cell-cell signaling. Required for normal granule cell migration in the developing cerebellum. Promotes reorganization of the actin cytoskeleton and plays an important role in axon guidance in the developing central nervous system.
Applications:	ELISA, WB
Name of antibody:	SEMA6A
Immunogen:	Synthetic peptide of humanSEMA6A

Full name:	sema domain, transmembrane domain (TM), and cytoplasmic domain, (semaphorin) 6A
Synonyms :	VIA; SEMA; HT018; SEMAQ; SEMA6A1
SwissProt:	Q9H2E6
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	114 kDa
WB Positive control:	293T, Lovo and HepG2 cells
WB Recommended dilution:	200-1000

