

兔抗 JAK2(Phospho-Tyr570) 多克隆抗体

- 中文名称： 兔抗 JAK2(Phospho-Tyr570) 多克隆抗体
- 英文名称： Anti-JAK2(Phospho-Tyr570) rabbit polyclonal antibody
- 别名： JTK10; THCYT3
- 相关类别： 一抗
- 储存： 冷冻（-20℃） 避光
- 宿主： Rabbit
- 抗原： JAK2(Phospho-Tyr570)
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.
Applications:	WB
Name of antibody:	JAK2(Phospho-Tyr570)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 570(G-D-Y(p)-G-Q) derived from Human JAK2 .

Full name:	Janus kinase 2
Synonyms :	JTK10; THCYT3
SwissProt:	O60674
WB Predicted band size :	131 kDa
WB Positive control:	Jurkat cells lysates
WB Recommended dilution:	500-1000

