

## 兔抗 VWA2 多克隆抗体

中文名称: 兔抗 VWA2 多克隆抗体

英文名称: Anti-VWA2 rabbit polyclonal antibody

别名: AMACO; NET42; CCSP-2

宿主: Rabbit

抗原: VWA2

标记物: Unconjugate

克隆类型: rabbit polyclonal

反应种属: Human

相关类别: 一抗 储存 冷冻 (-20℃)

### 技术规格

<b>Background:</b>	This gene encodes a member of the von Willebrand factor A-like domain protein superfamily. The encoded protein is localized to the extracellular matrix and may serve as a structural component in basement membranes or in anchoring structures on scaffolds of collagen VII or fibrillin. This gene has been linked to type 1A diabetes and is a candidate serological marker for colon cancer.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	VWA2
<b>Immunogen:</b>	Fusion protein of human VWA2
<b>Full name:</b>	von Willebrand factor A domain containing 2
<b>Synonyms:</b>	AMACO; NET42; CCSP-2

<b>SwissProt:</b>	Q5GFL6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	30-150

