

兔抗 SENP1 多克隆抗体

中文名称：兔抗 SENP1 多克隆抗体

英文名称：Anti-SENP1 rabbit polyclonal antibody

别名：SuPr-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SENP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a cysteine protease that specifically targets members of the small ubiquitin-like modifier (SUMO) protein family. This protease regulates SUMO pathways by deconjugating sumoylated proteins. This protease also functions to process the precursor SUMO proteins into their mature form. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	SENP1
Immunogen:	Synthetic peptide of human SENP1

Full name:	SUMO1/sentrin specific peptidase 1
Synonyms :	SuPr-2
SwissProt:	Q9P0U3
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	73 kDa
WB Positive control:	K562 cells
WB Recommended dilution:	200-1000



