

兔抗 SENP2 多克隆抗体

中文名称：兔抗 SENP2 多克隆抗体

英文名称： Anti-SENP2 rabbit polyclonal antibody

别名： AXAM2; SMT3IP2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SENP2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	SUMO1 (UBL1; MIM 601912) is a small ubiquitin-like protein that can be covalently conjugated to other proteins. SENP2 is one of a group of enzymes that process newly synthesized SUMO1 into the conjugatable form and catalyze the deconjugation of SUMO1-containing species.
Applications:	ELISA, IHC
Name of antibody:	SENP2
Immunogen:	Fusion protein of human SENP2
Full name:	SUMO1/sentrin/SMT3 specific peptidase 2
Synonyms :	AXAM2; SMT3IP2

SwissProt:	Q9HC62
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	25-100

