

兔抗 JAK2 多克隆抗体

中文名称: 兔抗 JAK2 多克隆抗体

英文名称: Anti-JAK2 rabbit polyclonal antibody

别名: JTK10, THCYT3

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: JAK2

反应种属: Human, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene product is a protein tyrosine kinase involved in a specific subset of cytokine receptor signaling pathways. It has been found to be constitutively associated with the prolactin receptor and is required for responses to gamma interferon. Mice that do not express an active protein for this gene exhibit embryonic lethality associated with the absence of definitive erythropoiesis.
Applications:	ELISA, IHC
Name of antibody:	JAK2
Immunogen:	Synthetic peptide of human JAK2

Full name:	Janus kinase 2
Synonyms :	JTK10, THCYT3
SwissProt:	O60674
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human lung cancer and Human ovarian cancer
IHC Recommend dilution:	10-50

