

兔抗 PDPK1 多克隆抗体

中文名称: 兔抗 PDPK1 多克隆抗体

英文名称: Anti-PDK1 (Ab-241) rabbit polyclonal antibody

别名: PDK1; PDPK2; PRO0461

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: PDPK1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Phosphorylates and activates not only PKB/AKT, but also PKA, PKC-zeta, RPS6KA1 and RPS6KB1. May play a general role in signaling processes and in development
Applications:	WB, IHC, IF
Name of antibody:	PDPK1
Immunogen:	Synthetic peptide of human PDK1 (Ab-241)
Full name:	3-phosphoinositide dependent protein kinase-1
Synonyms:	PDK1; PDPK2; PRO0461
SwissProt:	O15530
IHC positive control:	Human breast carcinoma

IHC Recommend dilution:	50-100
WB Predicted band size:	63 kDa
WB Positive control:	MCF, 293 and PC12 cells
WB Recommended dilution:	500-1000
IF Positive control:	MCF7 cells
IF Recommended dilution:	100-200



