

兔抗 VWF 多克隆抗体

中文名称：兔抗 VWF 多克隆抗体

英文名称： Anti-VWF rabbit polyclonal antibody

别名： VWD; F8VWF

储存： 冷冻（-20℃）

抗原： VWF

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The glycoprotein encoded by this gene functions as both an antihemophilic factor carrier and a platelet-vessel wall mediator in the blood coagulation system. It is crucial to the hemostasis process. Mutations in this gene or deficiencies in this protein result in von Willebrand's disease. An unprocessed pseudogene has been found on chromosome 22.
Applications:	ELISA, IHC
Name of antibody:	VWF

Immunogen:	Fusion protein of human VWF
Full name:	von Willebrand factor
Synonyms:	VWD; F8VWF
SwissProt:	P04275
IHC positive control:	Human lung cancer and human gastric cancer
IHC Recommend dilution:	25-100

