

## 兔抗 SENP7 多克隆抗体

中文名称：兔抗 SENP7 多克隆抗体

英文名称： Anti-SENP7 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SENP7

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The reversible posttranslational modification of proteins by the addition of small ubiquitin-like SUMO proteins (see SUMO1; MIM 601912) is required for many cellular processes. SUMO-specific proteases, such as SENP7, process SUMO precursors to generate a C-terminal diglycine motif required for the conjugation reaction.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SENP7
<b>Immunogen:</b>	Fusion protein of human SENP7
<b>Full name:</b>	SUMO1/sentrin specific peptidase 7
<b>SwissProt:</b>	Q9BQF6
<b>ELISA Recommended dilution:</b>	2000-10000

<b>WB Predicted band size:</b>	120 kDa
<b>WB Positive control:</b>	A172 and 231 cells
<b>WB Recommended dilution:</b>	1000-5000

