

兔抗 JAK3 多克隆抗体

中文名称: 兔抗 JAK3 多克隆抗体

英文名称: Anti-JAK3 rabbit polyclonal antibody

别名: JAK1; LAK; JAK-3; L-JAK; JAK3_HUMAN

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: JAK3

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the Janus kinase (JAK) family of tyrosine kinases involved in cytokine receptor-mediated intracellular signal transduction. It is predominantly expressed in immune cells and transduces a signal in response to its activation via tyrosine phosphorylation by interleukin receptors. Mutations in this gene are associated with autosomal SCID (severe combined immunodeficiency disease).
Applications:	ELISA, IHC
Name of antibody:	JAK3
Immunogen:	Synthetic peptide of human JAK3

Full name:	Janus kinase 3
Synonyms:	JAKL; LJAK; JAK-3; L-JAK; JAK3_HUMAN
SwissProt:	P52333
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

