

兔抗 JAKMIP2 多克隆抗体

中文名称： 兔抗 JAKMIP2 多克隆抗体

英文名称： Anti-JAKMIP2 rabbit polyclonal antibody

别名： NECC1; JAMIP2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： JAKMIP2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	JAKMIP2
Immunogen:	Fusion protein of human JAKMIP2
Full name:	janus kinase and microtubule interacting protein 2

Synonyms :	NECC1; JAMIP2
SwissProt:	Q96AA8
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human lung cancer and human brain
IHC Recommend dilution:	25-100

