

## 兔抗 SENP8 多克隆抗体

中文名称：兔抗 SENP8 多克隆抗体

英文名称： Anti-SENP8 rabbit ployclonal antibody

别名： DEN1; NEDP1; PRSC2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SENP8

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cysteine protease that is a member of the sentrin-specific protease family. The encoded protein is involved in processing and deconjugation of the ubiquitin-like protein termed, neural precursor cell expressed developmentally downregulated 8. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SENP8
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	SUMO/sentrin specific peptidase family member 8
<b>Synonyms :</b>	DEN1; NEDP1; PRSC2
<b>SwissProt:</b>	Q96LD8
<b>ELISA Recommended dilution:</b>	2000-5000

<b>IHC positive control:</b>	Human breast cancer and human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	24 kDa
<b>WB Positive control:</b>	293T, Hela, Lovo, Raji and huvec cells, Mouse spleen tissue
<b>WB Recommended dilution:</b>	500-2000



