

兔抗 SEPHS1 多克隆抗体

中文名称：兔抗 SEPHS1 多克隆抗体

英文名称： Anti-SEPHS1 rabbit polyclonal antibody

别名： selenophosphate synthetase 1; SPS; SELD; SPS1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEPHS1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an enzyme that synthesizes selenophosphate from selenide and ATP. Selenophosphate is the selenium donor used to synthesize selenocysteine, which is co-translationally incorporated into selenoproteins at in-frame UGA codons.
Applications:	ELISA, IHC
Name of antibody:	SEPHS1
Immunogen:	Fusion protein of human SEPHS1
Full name:	selenophosphate synthetase 1
Synonyms:	SPS; SELD; SPS1
SwissProt:	P49903

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	50-200

